DILATOR

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	7
Felhasználási javallatok	8
Indicazioni per l'uso	9
Bruksindikasjoner	10
Indicações de utilização	11
Indicaciones de uso	12
Indikationer för användning	13
Kullanım Endikasyonları	14

DILATOR

Indications for Use:

These dilators are used for the percutaneous introduction of guidewires into the peripheral vasculature.

Device Description:

The dilator introducer consists of a tipped radiopaque tube with a molded hub that accepts up to a .038 inch.

Contraindications:

Use of the dilator is contraindicated if the patient has a known or suspected obstruction in the vessel.

Potential Complications:

The potential complications related to the use of the dilator include, but are not limited to the following: Air embolism, device embolism, device dislodgement, pneumothorax, vein thrombosis, hematoma formation, hemothorax, vessel erosion, trauma to vessels, sepsis..

Precautions:

Store in a dry, dark, cool place. Do not use if package is open or damaged. Inspect all components prior to use.

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

Cautions:

- This procedure should only be performed by professionals thoroughly trained in this procedure.
- If resistance is met when advancing or withdrawing the guidewire through the needle, determine the cause by fluoroscopy and correct before continuing with the procedure.
- Do not attempt to use a guidewire over the maximum diameter specified on the package label.
- Individual patient anatomy and physician technique may require procedural variations.
- Insertion into artery may cause excessive bleeding and/or other complications.

Warnings:

- Do not alter this device in any way.
- Do not use alcohol, acetone or solutions containing these agents. These solutions may affect the properties of the plastic components resulting in degradation of the device.
- Do not reuse this device.
- Do not withdraw guidewire through metal needles; guidewire may shear or unravel.
- Do not resterilize.

USE STERILE TECHNIQUE. A suggested procedure:

1. Peel open package and place contents on sterile field.
2. Prep skin and drape in area of anticipated puncture site as desired.
3. Follow standard hospital practice for placement of dilator.

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Galt Medical, Corp.