Centeze®

Indications for Use	2
Indikationer for anvendelse	3
Gebruiksindicaties	4
Indications d'emploi	5
Anwendungsgebiete	6
Ενδείξεις χρήσης	7
Felhasználási javallatok	8
Indicazioni per l'uso	9
Bruksindikasjoner	10
Indicações de utilização	11
Indicaciones de uso	12
Indikationer för användning	13
Kullanım Endikasyonları	14

USA CAUTION

Federal (USA) law restricts this device to sale by or on the order of a physician.

DEVICE DESCRIPTION: The Centeze® centesis catheter (Centeze®) is an over-the needle catheter for percutaneous aspiration and small volume drainages for durations not to exceed 24 hours. Each Centeze® includes the following components: Catheter with three distal sideholes, percutaneous entry needle and porous plug.

INTENDED USE AND PATIENT POPULATION

The Centeze® is intended for use in percutaneous fluid aspirations and small volume drainage procedures. For use on both adult and adolescent groups.

CONTRAINDICATIONS

- Where percutaneous drainage catheterization is unacceptable
- Bleeding diathesis and uncontrolled hypertension
- Anticoagulant use
- Chemical ablation of tissue

WARNINGS

- The Centeze® is supplied sterile for single use only.
- Do not exceed insertion dwell time of 24 hours.
- Do not alter this device in any way.
- Do not resterilize the device. Re-sterilization could change the physical characteristics of the materials and should not be attempted.
- Do not reuse the device. Reuse will result in increased biocontamination risk for the patient resulting in infection or pyrogenic response.
- During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle and catheter tip.

PRECAUTIONS

The Centeze® deployment procedure should be performed by physicians thoroughly trained in drainage techniques and procedures.

Do not use the Centeze® if the packaging is damaged.

Inspect the Centeze® prior to use for any bends or kinks. Do not use a damaged needle because trauma to the patient and/or inability to place the drainage catheter may occur.

Any serious incident in relation to this device should be reported to the manufacturer and the competent authority established in the user and/or patient's Member State.

COMPLICATIONS: These are generally associated with drainage procedures and may occur when using the Centeze®:

- Local or systemic infection
- Pneumothorax
- Inadvertent arterial or venous puncture resulting in bleeding complications
- Catheter or needle fracture and embolization

CLINICAL PROCEDURE

The Centeze® should be used by physicians trained on the procedures for which the device is intended. The techniques and procedures described do not represent ALL medically acceptable protocols, nor are they intended as a substitute for the physician's experience and judgment in treating any specific patient. All available data, including the patient's signs and symptoms and other diagnostic test results, should be considered before determining a specific treatment plan.

PREPARATIONS FOR USE

1. Carefully inspect the Centeze® packaging and components for damage prior to use.
2. Using sterile technique, open the pouch and transfer the Centeze® into the sterile field.
3. Remove the needle from the catheter. Thoroughly flush both the needle and the catheter with saline, then reassemble them.
4. Prep skin and drape area of anticipated puncture as required.

DEPLOYMENT PROCEDURE

Deploy the Centeze® according to the following steps:

1. Prepare the skin incision site using local anesthetic.
2. Introduce the Centeze® into the cavity using ultrasound guidance, CT or fluoroscopy.
3. Remove the needle while leaving the catheter in position.
4. Connect a syringe or suction tube to the hub of the catheter to collect the drainage.
5. Upon completion of the drainage, remove the catheter and discard according to institutional practice.

WARNING: During insertion avoid contact with bone, cartilage and scar tissue, which can damage the needle and catheter tip.

PACKAGING & STORAGE

The Centeze® has been sterilized with ethylene oxide.

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