




GALT™

OEM Division

A Theragenics Company



Galt Vascular and Interventional Medical Devices

Product Catalog

www.galtmedical.com

Welcome

Galt Medical Corp. has been a leading provider of vascular and interventional medical devices for over 30 years, developing and manufacturing a range of products for interventional procedures. Our regulatory, customer service, engineering, and marketing teams are located in our Garland, Texas headquarters, with manufacturing done both at our Garland location and in Costa Rica.

Galt products are widely recognized by our OEM customers and clinicians for their high-quality design, consistency, reliability, and value. Our manufacturing capabilities include, but are not limited to tipping, forming, shaping, injection molding, hot stamping, hole punching, markings, assembly, packaging, tooling, and rapid prototyping/3D printing. Galt also provides regulatory support to many of its contract manufacturing customers.

Projects range from exclusive, minor modifications of existing Galt products for an OEM partner, to unique products based on an OEM's requirements or specifications. We listen and work closely with our end user and OEM customers to make sure our products address real patient and partner needs. While many of our products serve broad markets and utilize our high-volume manufacturing capabilities, we are proud that several of our products address under-served, niche applications.

Get to know Galt Medical Corp. and our complete offering of access devices. From full kits to individual components, our products are available bulk non-sterile or Galt private labeled.

Trusted Experience. Industry Leaders.

Galt Medical Corp. and TGXCR (Costa Rica) is FDA registered, ISO 13485:2016, EN/ISO 13485:2016, Directive 93/42/EEC certified. We have MDSAP certifications in Australia, Brazil, Canada, and United States. MDD class II and III certifications. Currently our MDR certification is pending. Please contact us for updates.



Contents

Guidewires & Accessories



Page 4

Valved Tearaway Introducers



Page 13

Hemostasis Valve Introducers & Accessories




Page 8

Micro Introducer Kits & Sets



Page 14

Tearaway Introducers




Page 11

Needles



Page 16



Guidewires

SKU	Product Description	Outer Diameter	Length	Core Material	Coil Material	Safety Wire	Tip Configuration	UoM
GWA-015-21	.018 x 40cm, SS/SS	0.018	40cm	Stainless Steel	Stainless Steel	No	Straight	1 ea
GWA-015-02	.018 x 40cm, SS/SS	0.018	40cm	Stainless Steel	Stainless Steel	No	Straight	1 ea
GWA-055-41	.018 x 40cm, SS/AU	0.018	40cm	Stainless Steel	Gold	No	Straight	1 ea
GWA-055-07	.018 x 40cm, SS/AU	0.018	40cm	Stainless Steel	Gold	No	Straight	1 ea
GWA-015-07	.018 x 40cm, SS/PT	0.018	40cm	Stainless Steel	Platinum	No	Straight	1 ea
GWA-019-30	.018 x 40cm, NT/SS	0.018	40cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-019-52	.018 x 40cm, NT/SS	0.018	40cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-019-29	.018 x 40cm, NT/SS, Angled Tip	0.018	40cm	Nitinol	Stainless Steel	No	Angled	1 ea
GWA-056-46	.018 x 40cm, NT/AU	0.018	40cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-01	.018 x 40cm, NT/AU	0.018	40cm	Nitinol	Gold	No	Straight	1 ea
GWA-019-10	.018 x 40cm, NT/PT	0.018	40cm	Nitinol	Platinum	No	Straight	1 ea
GWA-019-31	.018 x 45cm, NT/SS	0.018	45cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-022-69	.018 x 45cm, NT/SS	0.018	45cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-056-04	.018 x 45cm, NT/AU	0.018	45cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-11	.018 x 45cm, NT/AU, Angled Tip	0.018	45cm	Nitinol	Gold	No	Angled	1 ea
GWA-022-37	.018 x 45cm, NT/PT	0.018	45cm	Stainless Steel	Platinum	No	Straight	1 ea
GWA-022-67	.018 x 50cm, NT/SS, Marked	0.018	50cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-038-24	.018 x 50cm, NT/AU, Angled Tip	0.018	80cm	Nitinol	Gold	No	Angled - 15 °	1 ea
GWA-015-03	.018 x 60cm, SS/SS	0.018	60cm	Stainless Steel	Stainless Steel	No	Straight	1 ea
GWA-055-42	.018 x 60cm, SS/AU	0.018	60cm	Stainless Steel	Gold	No	Straight	1 ea
GWA-019-32	.018 x 60cm, NT/SS	0.018	60cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-056-48	.018 x 60cm, NT/AU	0.018	60cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-02	.018 x 60cm, NT/AU	0.018	60cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-06	.018 x 60cm, NT/AU	0.018	60cm	Nitinol	Gold	No	Straight	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label

Guidewires continued...

SKU	Product Description	Outer Diameter	Length	Core Material	Coil Material	Safety Wire	Tip Configuration	UoM
GWA-022-41	.018 x 65cm, SS/PT	0.018	65cm	Stainless Steel	Platinum	No	Straight	1 ea
GWA-022-11	.018 x 65cm, NT/SS, Angled Tip	0.018	65cm	Nitinol	Stainless Steel	No	Angled	1 ea
GWA-019-33	.018 x 70cm, NT/SS	0.018	70cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-056-49	.018 x 70cm, NT/AU	0.018	70cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-38	.018 x 70cm, NT/AU, Marked	0.018	70cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-39	.018 x 70cm, NT/AU, Marked	0.018	70cm	Nitinol	Gold	No	Straight	1 ea
GWA-015-14	.018 x 80cm, SS/SS	0.018	80cm	Stainless Steel	Stainless Steel	No	Straight	1 ea
GWA-056-03	.018 x 80cm, NT/AU	0.018	80cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-53	.018 x 80cm, NT/AU	0.018	80cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-43	.018 x 80cm, NT/AU, Marked	0.018	80cm	Nitinol	Gold	No	Straight	1 ea
GWA-019-14	.018 x 80cm, NT/PT	0.018	80cm	Nitinol	Platinum	No	Straight	1 ea
GWA-019-26	.018 x 80cm, NT/PT, Angled Tip	0.018	80cm	Nitinol	Platinum	No	Angled	1 ea
GWA-019-34	.018 x 80cm, NT/SS	0.018	80cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-019-42	.018 x 100cm, NT/SS	0.018	100cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-019-43	.018 x 120cm, NT/SS	0.018	120cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-022-09	.018 x 130cm, SS/SS	0.018	130cm	Stainless Steel	Stainless Steel	No	Straight	1 ea
GWA-022-13	.018 x 130cm, NT/SS	0.018	130cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-056-40	.018 x 130cm, NT/AU, Marked	0.018	130cm	Nitinol	Gold	No	Straight	1 ea
GWA-056-41	.018 x 130cm, NT/AU, Marked	0.018	130cm	Nitinol	Gold	No	Straight	1 ea
GWA-019-44	.018 x 150cm, NT/SS	0.018	150cm	Nitinol	Stainless Steel	No	Straight	1 ea
GWA-038-25	.018 x 180cm, NT/AU, TORX	0.018	180cm	Nitinol	Gold	No	Straight	1 ea
GWA-038-20	.018 x 300cm, NT/AU, TORX	0.018	300cm	Nitinol	Gold	No	Straight	1 ea
GWA-025-06	.035 x 50cm, SS/SS	0.0350	50cm	Stainless Steel	Stainless Steel	Yes	3.0 mm J Tip	1 ea
GWA-012-09	.035 x 100cm, SS/SS	0.035	100cm	Stainless Steel	Stainless Steel	Yes	3.0 mm J Tip	1 ea
GWA-001-15	.035 x 150cm, SS/SS	0.035	150cm	Stainless Steel	Stainless Steel	Yes	Straight	1 ea
GWA-012-10	.035 x 150cm, SS/SS	0.035	150cm	Stainless Steel	Stainless Steel	Yes	3.0 mm J Tip	1 ea
GWA-004-13	.038 x 50cm, SS/SS	0.038	50cm	Stainless Steel	Stainless Steel	Yes	3.0 mm J Tip	1 ea
GWA-003-13	.038 x 70cm, SS/SS	0.038	70cm	Stainless Steel	Stainless Steel	Yes	3.0 mm J Tip	1 ea
GWA-012-16	.038 x 100cm, SS/SS	0.038	100cm	Stainless Steel	Stainless Steel	Yes	3.0 mm J Tip	1 ea
GWA-001-25	.038 x 150cm, SS/SS	0.038	150cm	Stainless Steel	Stainless Steel	Yes	Straight	1 ea
GWA-038-25	.018 x 180cm, NT/AU, TORX	0.018	180cm	Nitinol	Gold	No	Straight	1 ea
GWA-038-20	.018 x 300cm, NT/AU, TORX	0.018	300cm	Nitinol	Gold	No	Straight	1 ea
GWA-019-17	.018 x 300cm, NT/PT, Angled Tip	0.018	300cm	Nitinol	Platinum	No	Angled	1 ea
GWA-004-07	.035 x 70cm, SS/SS	0.035	70cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-004-08	.035 x 80cm, SS/SS	0.035	80cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-012-09	.035 x 100cm, SS/SS	0.035	100cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-001-15	.035 x 150cm, SS/SS	0.035	150cm	Stainless Steel	Stainless Steel	Yes	Straight	1 ea
GWA-012-10	.035 x 150cm, SS/SS	0.035	150cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-004-13	.038 x 50cm, SS/SS	0.038	50cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-003-13	.038 x 70cm, SS/SS	0.038	70cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-012-16	.038 x 100cm, SS/SS	0.038	100cm	Stainless Steel	Stainless Steel	Yes	3.0cm J Tip	1 ea
GWA-001-25	.038 x 150cm, SS/SS	0.038	150cm	Stainless Steel	Stainless Steel	Yes	Straight	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label

Guidewires in Dispenser


SKU	Product Description	Outer Diameter	Length	Core Material	Coil Material	Tip Configuration	UoM
GID-015-21	.018 x 40cm, SS/SS	0.018	40cm	Stainless Steel	Stainless Steel	Straight	1 ea
GID-015-07	.018 x 40cm, SS/PT	0.018	40cm	Stainless Steel	Platinum	Straight	1 ea
GID-019-30	.018 x 40cm, NT/SS	0.018	40cm	Nitinol	Stainless Steel	Straight	1 ea
GID-056-46	.018 x 40cm, NT/AU	0.018	40cm	Nitinol	Gold	Straight	1 ea
GID-019-10	.018 x 40cm, NT/PT	0.018	40cm	Nitinol	Platinum	Straight	1 ea
GID-019-31	.018 x 45cm, NT/SS	0.018	45cm	Nitinol	Stainless Steel	Straight	1 ea
GID-056-04	.018 x 45cm, NT/AU	0.018	45cm	Nitinol	Gold	Straight	1 ea
GID-056-11	.018 x 45cm, NT/AU, Angled Tip	0.018	45cm	Nitinol	Gold	Angled	1 ea
GID-022-67	.018 x 50cm, NT/SS, Marked	0.018	50cm	Nitinol	Stainless Steel	Straight	1 ea
GID-055-42	.018 x 60cm, SS/AU	0.018	60cm	Stainless Steel	Gold	Straight	1 ea
GID-019-32	.018 x 60cm, NT/SS	0.018	60cm	Nitinol	Stainless Steel	Straight	1 ea
GID-056-06	.018 x 60cm, NT/AU	0.018	60cm	Nitinol	Gold	Straight	1 ea
GID-022-41	.018 x 65cm, SS/PT, Marked	0.018	65cm	Stainless Steel	Platinum	Straight	1 ea
GID-022-11	.018 x 65cm, NT/SS, Angled Tip	0.018	65cm	Nitinol	Stainless Steel	Angled	1 ea
GID-019-33	.018 x 70cm, NT/SS	0.018	70cm	Nitinol	Stainless Steel	Straight	1 ea
GID-056-38	.018 x 70cm, NT/AU, Marked	0.018	70cm	Nitinol	Gold	Straight	1 ea
GID-056-39	.018 x 70cm, NT/AU, Marked	0.018	70cm	Nitinol	Gold	Straight	1 ea
GID-019-34	.018 x 80cm, NT/SS	0.018	80cm	Nitinol	Stainless Steel	Straight	1 ea
GID-056-43	.018 x 80cm, NT/AU, Marked	0.018	80cm	Nitinol	Gold	Straight	1 ea
GID-019-14	.018 x 80cm, NT/PT	0.018	80cm	Nitinol	Platinum	Straight	1 ea
GID-019-26	.018 x 80cm, NT/PT, Angled Tip	0.018	80cm	Nitinol	Platinum	Angled	1 ea
GID-019-42	.018 x 100cm, NT/SS	0.018	100cm	Nitinol	Stainless Steel	Straight	1 ea
GID-019-43	.018 x 120cm, NT/SS	0.018	120cm	Nitinol	Stainless Steel	Straight	1 ea
GID-022-13	.018 x 130cm, NT/SS, Marked	0.018	130cm	Nitinol	Stainless Steel	Straight	1 ea
GID-056-40	.018 x 130cm, NT/AU, Marked	0.018	130cm	Nitinol	Gold	Straight	1 ea
GID-056-41	.018 x 130cm, NT/AU, Marked	0.018	130cm	Nitinol	Gold	Straight	1 ea
GID-012-09	.035 x 100cm, SS/SS	0.035	100cm	Stainless Steel	Stainless Steel	3.0 mm J Tip	1 ea
GID-012-10	.035 x 150cm, SS/SS	0.035	150cm	Stainless Steel	Stainless Steel	3.0 mm J Tip	1 ea
GID-004-13-01	.038 x 50cm, SS/SS	0.038	50cm	Stainless Steel	Stainless Steel	3.0 mm J Tip	1 ea
GID-003-13S	.038 x 70cm, SS/SS	0.038	70cm	Stainless Steel	Stainless Steel	3.0 mm J Tip	1 ea
GID-012-16	.038 x 100cm, SS/SS	0.038	100cm	Stainless Steel	Stainless Steel	3.0 mm J Tip	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label



Twisted Wire Stylet

SKU	Product Description	Outer Diameter	Length	Core Material	Coil Material	Safety Wire	Tip Plating	Tip Configuration	UoM
STY-000-07	.009 x 72cm Twisted Wire Stylet	0.009	72cm	n/a	n/a	n/a	n/a	n/a	1 ea
STY-000-24	.012 x 75cm Twisted Wire Stylet	0.012	75cm	n/a	n/a	n/a	n/a	n/a	1 ea
STY-001-25	.012 x 85cm Twisted Wire Stylet	0.012	85cm	n/a	n/a	n/a	n/a	n/a	1 ea
STY-000-37	.015 x 37cm Twisted Wire Stylet	0.015	37cm	n/a	n/a	n/a	n/a	n/a	1 ea
STY-000-43	.015 x 80cm Twisted Wire Stylet	0.015	80cm	n/a	n/a	n/a	n/a	n/a	1 ea



Torque Device

SKU	Product Description	Outer Diameter	Length	Core Material	Coil Material	Safety Wire	Tip Plating	Tip Configuration	UoM
TOR-000-10	.010 - .020 QuikTork (Yellow)	n/a	n/a	n/a	n/a	n/a	n/a	n/a	1 ea
TOR-000-25	.025 - .040 QuikTork (Orange)	n/a	n/a	n/a	n/a	n/a	n/a	n/a	1 ea

Hemostasis Valve Introducers & Accessories




Hemostasis Valve Introducers

SKU	Product Description	Sheath Diameter	Sheath Length	Dilator Length	Dilator Tipped	Coating	UoM
CLI-200-21	Cath Lab Sheath with Side Port	4.0F	12cm	16.25cm	0.019		1 ea
CLI-200-30	Cath Lab Sheath with Side Port	4.0F	24cm	28.25cm	0.019		1 ea
CLI-200-31	Cath Lab Sheath with Side Port	4.0F	45cm	49.25cm	0.019		1 ea
CLI-200-32	Cath Lab Sheath with Side Port	4.0F	55cm	59.25cm	0.019		1 ea
CLI-200-22	Cath Lab Sheath with Side Port	5.5F	12cm	16.25cm	0.019		1 ea
CLI-200-35	Cath Lab Sheath with Side Port	5.5F	24cm	30.30cm	0.019		1 ea
CLI-200-26	Cath Lab Sheath with Side Port	6.5F	12cm	18.30cm	0.019		1 ea
CLI-200-36	Cath Lab Sheath with Side Port	6.5F	24cm	30.30cm	0.019		1 ea
CLI-200-37	Cath Lab Sheath with Side Port	7.0F	24cm	30.30cm	0.019		1 ea
CLI-202-08	Cath Lab Sheath with Side Port	8.0F	12cm	17.50cm	0.039		1 ea
CLI-202-18	Cath Lab Sheath with Side Port	8.0F	24cm	29.50cm	0.039		1 ea
CLI-210-34	Cath Lab Sheath with Side Port	4.0F	10cm	14.25cm	0.019	Silicone	1 ea
CLI-210-36	Cath Lab Sheath with Side Port	4.0F	12cm	16.25cm	0.019	Silicone	1 ea
CLI-210-37	Cath Lab Sheath with Side Port	4.0F	24cm	28.25cm	0.019	Silicone	1 ea
CLI-210-15	Cath Lab Sheath with Side Port	5.5F	7cm	11.05cm	0.019	Silicone	1 ea
CLI-212-05	Cath Lab Sheath with Side Port	5.5F	12cm	16.25cm	0.039	Silicone	1 ea
CLI-212-15	Cath Lab Sheath with Side Port	5.5F	24cm	29.50cm	0.039	Silicone	1 ea
CLI-210-23	Cath Lab Sheath with Side Port	6.0F	12cm	16.25cm	0.019	Silicone	1 ea
CLI-210-16	Cath Lab Sheath with Side Port	6.5F	7cm	11.05cm	0.019	Silicone	1 ea
CLI-212-06	Cath Lab Sheath with Side Port	6.5F	12cm	16.25cm	0.039	Silicone	1 ea
CLI-212-16	Cath Lab Sheath with Side Port	6.5F	24cm	29.50cm	0.039	Silicone	1 ea
CLI-210-27	Cath Lab Sheath with Side Port	7.0F	7cm	13.30cm	0.019	Silicone	1 ea
CLI-210-28	Cath Lab Sheath with Side Port	7.0F	12cm	16.25cm	0.019	Silicone	1 ea
CLI-212-17	Cath Lab Sheath with Side Port	7.0F	24cm	29.50cm	0.039	Silicone	1 ea
CLI-210-17	Cath Lab Sheath with Side Port	8.0F	7cm	11.25cm	0.019	Silicone	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label | Sheaths and dilators can be sold separately

Hemostasis Valve Introducers continued...

SKU	Product Description	Sheath Diameter	Sheath Length	Dilator Length	Dilator Tipped	Coating	UoM
CLI-210-29	Cath Lab Sheath with Side Port	8.0F	12cm	16.25cm	0.019	Silicone	1 ea
CLI-210-35	Cath Lab Sheath with Side Port	8.0F	24cm	28.25cm	0.019	Silicone	1 ea
CLI-212-18	Cath Lab Sheath with Side Port	8.0F	24cm	29.50cm	0.039	Silicone	1 ea
CLI-212-09	Cath Lab Sheath with Side Port	9.0F	12cm	17.50cm	0.039	Silicone	1 ea



Radial Hemostasis Valve Introducers

SKU	Product Description	Sheath Diameter	Sheath Length	Dilator Length	Dilator Tipped	Coating	Type of Wire	Type of Needle	UoM
CLI-228-45	Hydrophilic Radial Artery Access Kit	4.5F	7cm	11.8cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-230-45	Hydrophilic Radial Artery Access Kit	4.5F	7cm	11.8cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-229-45	Hydrophilic Radial Artery Access Kit	4.5F	11cm	16.3cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-231-45	Hydrophilic Radial Artery Access Kit	4.5F	11cm	16.3cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-228-55	Hydrophilic Radial Artery Access Kit	5.5F	7cm	11.8cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-230-55	Hydrophilic Radial Artery Access Kit	5.5F	7cm	11.8cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-229-55	Hydrophilic Radial Artery Access Kit	5.5F	11cm	16.6cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-231-55	Hydrophilic Radial Artery Access Kit	5.5F	11cm	16.6cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-228-06	Hydrophilic Radial Artery Access Kit	6.0F	7cm	12.3cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-230-06	Hydrophilic Radial Artery Access Kit	6.0F	7cm	12.3cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-229-06	Hydrophilic Radial Artery Access Kit	6.0F	11cm	15.8cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box

Radial Artery Kits: Can only be sold as sterile with your own company private label

Hemostasis Valve Introducers & Accessories

Radial Hemostasis Valve Introducers continued...

SKU	Product Description	Sheath Diameter	Sheath Length	Dilator Length	Dilator Tipped	Coating	Type of Wire	Type of Needle	UoM
CLI-231-06	Hydrophilic Radial Artery Access Kit	6.0F	11cm	15.8cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-228-07	Hydrophilic Radial Artery Access Kit	7.0F	7cm	12.6cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-230-07	Hydrophilic Radial Artery Access Kit	7.0F	7cm	12.6cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-229-07	Hydrophilic Radial Artery Access Kit	7.0F	11cm	15.8cm	0.019	Hydrophilic	.018 x 45cm, 7cm SS Coil	21g x 2.5cm Introducer Needle, Non Echo	10ea / box
CLI-231-07	Hydrophilic Radial Artery Access Kit	7.0F	11cm	15.8cm	0.019	Hydrophilic	.018 x 45cm, NT	21g x 2.5cm Introducer Needle, Non Echo	10ea / box



Radial Hemostasis Valve Adaptor

Hemostasis Valve Adaptor

SKU	Product Description	Sheath Diameter	Sheath Length	Dilator Length	Dilator Tipped	Coating	UoM
ACL-001-00	Hemostasis Valve Adaptor	n/a	n/a	n/a	n/a	n/a	1 ea

Can be sold in either BNS and/or sterile private label



Tearaway Sheath / Dilator

SKU	Product Description	Sheath Diameter	Sheath Length	Sheath Material	Dilator Length	Dilator Tipped	UoM
DSS-005-05	T/A Int Set 5F x 14cm	5.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-06	T/A Int Set 6F x 14cm	6.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-065	T/A Int Set 6.5F x 14cm	6.5F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-012-07	T/A Int Set 7F x 12cm	7.0F	12cm	PTFE	16.7cm	0.039	1 ea
DSS-005-07	T/A Int Set 7F x 14cm	7.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-08	T/A Int Set 8F x 14cm	8.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-012-09	T/A Int Set 9F x 12cm	9.0F	12cm	PTFE	16.7cm	0.039	1 ea
DSS-005-09	T/A Int Set 9F x 14cm	9.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-020-09	T/A Int Set 9F x 20cm	9.0F	20cm	PTFE	23.6cm	0.039	1 ea
DSS-005-10	T/A Int Set 10F x 14cm	10.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-105	T/A Int Set 10.5F x 14cm	10.5F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-11	T/A Int Set 11F x 14cm	11.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-12	T/A Int Set 12F x 14cm	12.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-13	T/A Int Set 13F x 14cm	13.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-005-14	T/A Int Set 14F x 14cm	14.0F	14cm	PTFE	17.75cm	0.039	1 ea
DSS-020-15	T/A Int Set 15F x 20cm	15.0F	20cm	PTFE	23.6cm	0.039	1 ea
DSS-005-16	T/A Int Set 16F x 14cm	16.0F	14cm	PTFE	17.75cm	0.039	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label | Sheaths and dilators can be sold separately

Tearway Introducers



Micro Access Tearway Sheath / Dilator

SKU	Product Description	Sheath Diameter	Sheath Length	Sheath Material	Dilator Length	Dilator Tipped	UoM
DSS-024-15	Micro Access 2F x 7cm	2.0F	7cm	PTFE	8.80cm	0.019	1 ea
DSS-024-16	Micro Access 2F x 10cm	2.0F	10cm	PTFE	11.80cm	0.019	1 ea
DSS-006-03	Micro Access 3F x 5cm	3.0F	5cm	PTFE	7.50cm	0.019	1 ea
DSS-010-03	Micro Access 3F x 7cm	3.0F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-03	Micro Access 3F x 10cm	3.0F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-006-035	Micro Access 3.5F x 5cm	3.5F	5cm	PTFE	7.50cm	0.019	1 ea
DSS-010-035	Micro Access 3.5F x 7cm	3.5F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-035	Micro Access 3.5F x 10cm	3.5F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-006-04	Micro Access 4F x 5cm	4.0F	5cm	PTFE	7.50cm	0.019	1 ea
DSS-010-04	Micro Access 4F x 7cm	4.0F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-04	Micro Access 4F x 10cm	4.0F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-006-045	Micro Access 4.5F x 5cm	4.5F	5cm	PTFE	7.50cm	0.019	1 ea
DSS-010-045	Micro Access 4.5F x 7cm	4.5F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-045	Micro Access 4.5F x 10cm	4.5F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-006-05	Micro Access 5F x 5cm	5.0F	5cm	PTFE	7.95cm	0.019	1 ea
DSS-010-05	Micro Access 5F x 7cm	5.0F	7cm	PTFE	9.70cm	0.019	1 ea
DSS-007-05	Micro Access 5F x 10cm	5.0F	10cm	PTFE	12.70cm	0.019	1 ea
DSS-006-055	Micro Access 5.5F x 5cm	5.5F	5cm	PTFE	7.50cm	0.019	1 ea
DSS-010-055	Micro Access 5.5F x 7cm	5.5F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-055	Micro Access 5.5F x 10cm	5.5F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-006-06	Micro Access 6F x 5cm	6.0F	5cm	PTFE	7.50cm	0.019	1 ea
DSS-010-06	Micro Access 6F x 7cm	6.0F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-06	Micro Access 6F x 10cm	6.0F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-010-065	Micro Access 6.5F x 7cm	6.5F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-065	Micro Access 6.5F x 10cm	6.5F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-010-07	Micro Access 7F x 7cm	7.0F	7cm	PTFE	9.25cm	0.019	1 ea
DSS-007-07	Micro Access 7F x 10cm	7.0F	10cm	PTFE	12.25cm	0.019	1 ea
DSS-010-10	Micro Access 10F x 7cm	10.0F	7cm	PTFE	9.25cm	0.019	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label | Sheaths and dilators can be sold separately




Valved Tearaway Sheath / Dilator

SKU	Product Description	Sheath Diameter	Sheath Length	Sheath Material	Dilator Length	Dilator Tipped	Additional Features	UoM
INT-042-05	Valved T/A Int Set 5F x 13cm	5.0F	13cm	PTFE	16.60cm	0.039		1 ea
INT-042-06	Valved T/A Int Set 6F x 13cm	6.0F	13cm	PTFE	16.60cm	0.039		1 ea
INT-042-07	Valved T/A Int Set 7F x 13cm	7.0F	13cm	PTFE	16.70cm	0.039		1 ea
INT-042-08	Valved T/A Int Set 8F x 13cm	8.0F	13cm	PTFE	17.00cm	0.039		1 ea
INT-042-09	Valved T/A Int Set 9F x 13cm	9.0F	13cm	PTFE	17.20cm	0.039		1 ea
INT-042-10	Valved T/A Int Set 10F x 13cm	10.0F	13cm	PTFE	17.60cm	0.039		1 ea
INT-042-11	Valved T/A Int Set 11F x 13cm	11.0F	13cm	PTFE	17.60cm	0.039		1 ea
INT-042-12	Valved T/A Int Set 12F x 13cm	12.0F	13cm	PTFE	17.60cm	0.039		1 ea
INT-042-13	Valved T/A Int Set 13F x 13cm	13.0F	13cm	PTFE	17.60cm	0.039		1 ea
INT-042-14	Valved T/A Int Set 14F x 13cm	14.0F	13cm	PTFE	17.60cm	0.039		1 ea
INT-042-15	Valved T/A Int Set 15F x 13cm	15.0F	13cm	PTFE	17.60cm	0.039		1 ea
INT-042-16	Valved T/A Int Set 16F x 13cm	16.0F	13cm	PTFE	17.70cm	0.039		1 ea
INT-043-05	Valved T/A Int Set with Side Arm/Port 5F x 13cm	5.0F	13cm	PTFE	16.60cm	0.039	Side Arm/Port	1 ea
INT-043-06	Valved T/A Int Set with Side Arm/Port 6F x 13cm	6.0F	13cm	PTFE	16.60cm	0.039	Side Arm/Port	1 ea
INT-043-07	Valved T/A Int Set with Side Arm/Port 7F x 13cm	7.0F	13cm	PTFE	16.70cm	0.039	Side Arm/Port	1 ea
INT-043-08	Valved T/A Int Set with Side Arm/Port 8F x 13cm	8.0F	13cm	PTFE	17.00cm	0.039	Side Arm/Port	1 ea
INT-043-09	Valved T/A Int Set with Side Arm/Port 9F x 13cm	9.0F	13cm	PTFE	17.20cm	0.039	Side Arm/Port	1 ea
INT-043-10	Valved T/A Int Set with Side Arm/Port 10F x 13cm	10.0F	13cm	PTFE	17.60cm	0.039	Side Arm/Port	1 ea
INT-043-11	Valved T/A Int Set with Side Arm/Port 11F x 13cm	11.0F	13cm	PTFE	17.60cm	0.039	Side Arm/Port	1 ea
INT-043-12	Valved T/A Int Set with Side Arm/Port 12F x 13cm	12.0F	13cm	PTFE	17.60cm	0.039	Side Arm/Port	1 ea
INT-043-13	Valved T/A Int Set with Side Arm/Port 13F x 13cm	13.0F	13cm	PTFE	17.60cm	0.039	Side Arm/Port	1 ea
INT-043-14	Valved T/A Int Set with Side Arm/Port 14F x 13cm	14.0F	13cm	PTFE	17.60cm	0.039	Side Arm/Port	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label | Sheaths and dilators can be sold separately

Micro Introducer Kits & Sets



Micro Introducer Kits

SKU	Product Description	Sheath Diameter	Sheath Length	Type of Wire	Type of Needle	UoM
KIT-002-30	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, SS/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-26	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, SS/AU	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-28	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, SS/AU	21g x 7cm, Introducer Needle, Non Echo	10ea / box
KIT-002-32	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, NT/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-34	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, NT/AU	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-36	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 60cm, SS/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-38	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 60cm, NT/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-40	Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 60cm, NT/AU	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-31	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, SS/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-27	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, SS/AU	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-29	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, SS/AU	21g x 7cm, Introducer Needle, Non Echo	10ea / box
KIT-002-33	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, NT/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-35	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, NT/AU	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-37	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 60cm, SS/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-002-39	Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 60cm, NT/SS	21g x 7cm, Introducer Needle, Echo	10ea / box

Stiffened Micro Introducer Kits

SKU	Product Description	Sheath Diameter	Sheath Length	Type of Wire	Type of Needle	UoM
KIT-011-02	Stiffened Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, SS/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-011-20	Stiffened Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, SS/AU	21g x 7cm, Introducer Needle, Non Echo	10ea / box
KIT-011-40	Stiffened Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, NT/SS	21g x 7cm, Introducer Needle, Non Echo	10ea / box
KIT-011-42	Stiffened Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, NT/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-011-62	Stiffened Micro-Introducer Kit, 4F x 10cm	4.0F	10cm	.018 x 40cm, NT/AU	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-011-99	Stiffened Micro-Introducer Kit, 4F x 15cm	4.0F	15cm	.018 x 60cm, NT/AU	21g x 12cm, Introducer Needle, Echo	10ea / box
KIT-011-03	Stiffened Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, SS/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-011-41	Stiffened Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, NT/SS	21g x 7cm, Introducer Needle, Non Echo	10ea / box
KIT-011-43	Stiffened Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, NT/SS	21g x 7cm, Introducer Needle, Echo	10ea / box
KIT-011-63	Stiffened Micro-Introducer Kit, 5F x 10cm	5.0F	10cm	.018 x 40cm, NT/AU	21g x 7cm, Introducer Needle, Echo	10ea / box

Micro Introducer Kits : Custom sizes and modifications available | Can only be sold as sterile with your own company private label



GALTstick						
SKU	Product Description	Sheath Diameter	Sheath Length	Type of Wire	Type of Needle	UoM
KIT-030-06	6F GALTstick Kit	6.0F	19.8cm	.018 x 60cm, SS/AU	21g x 15cm, Trocar	10ea / box
KIT-030-03	6F GALTstick Kit	6.0F	19.8cm	.018 x 60cm, NT/AU	21g x 15cm, Trocar	10ea / box
KIT-036-06	6F Hydrophilic GALTstick Kit	6.0F	19.8cm	.018 x 60cm, SS/AU	21g x 15cm, Trocar	10ea / box
KIT-036-03	6F Hydrophilic GALTstick Kit	6.0F	19.8cm	.018 x 60cm, NT/AU	21g x 15cm, Trocar	10ea / box

GALTstick : Can only be sold as sterile with your own company private label



Introducer Needles

SKU	Product Description	Gauge	Length	Non Echo/Echo	Bevel	Hub Color	Wall Diameter Description	UoM
NDL-000-00	18g x 7cm, Introducer Needle	18g	7cm	Non Echo	B Bevel	Pink Hub	Extra Thin Wall	1 ea
NDL-000-33	18g x 7cm, Introducer Needle	18g	7cm	Echo	B Bevel	Pink Hub	Extra Thin Wall	1 ea
NDL-000-02	18g x 9cm, Introducer Needle	18g	9cm	Echo	B Bevel	Pink Hub	Extra Thin Wall	1 ea
NDL-000-07	19g x 5cm, Introducer Needle	19g	5cm	Non Echo	B Bevel	Brown Hub	Extra Thin Wall	1 ea
NDL-000-09	19g x 7cm, Introducer Needle	19g	7cm	Non Echo	B Bevel	Brown Hub	Extra Thin Wall	1 ea
NDL-000-10	19g x 7cm, Introducer Needle	19g	7cm	Echo	B Bevel	Brown Hub	Extra Thin Wall	1 ea
NDL-000-15	21g x 2cm Introducer Needle	21G	2cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-32	21g x 2.5cm, Introducer Needle	21g	2.5cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-17	21g x 4cm, Introducer Needle	21g	4cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-007-14	21g x 4cm, X-Sharp Introducer Needle	21g	4cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-18	21g x 4cm, Introducer Needle	21g	4cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-007-04	21g x 4cm, X-Sharp Introducer Needle	21g	4cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-20	21g x 5cm, Introducer Needle	21g	5cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-21	21g x 7cm, Introducer Needle	21g	7cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-007-17	21g x 7cm, X-Sharp Introducer Needle	21g	7cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-22	21g x 7cm, Introducer Needle	21g	7cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-007-07	21g x 7cm, X-Sharp Introducer Needle	21g	7cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-23	21g x 9cm, Introducer Needle	21g	9cm	Non Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-24	21g x 9cm, Introducer Needle	21g	9cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-28	21g x 12cm, Introducer Needle, Echo	21G	12cm	Echo	B Bevel	Green Hub	Thin Wall	1 ea
NDL-000-29	24g x 2.5cm, Introducer Needle	24g	2.5cm	Echo	B Bevel	Clear Hub	Extra Thin Wall	1 ea
NDL-000-30	24g x 4cm, Introducer Needle	24g	4cm	Echo	B Bevel	Clear Hub	Extra Thin Wall	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label



Centesis Catheter Needles

SKU	Product Description	Gauge	Length	Non Echo/Echo	Bevel	Hub Color	Wall Diameter Description	UoM
NDL-009-04	17g x 9cm, Centesis Catheter Needle	17g	9cm	n/a	n/a	Clear Hub	n/a	1 ea
NDL-009-00	17g x 11cm, Centesis Catheter Needle	17g	11cm	n/a	n/a	Clear Hub	n/a	1 ea
NDL-009-05	17g x 12cm, Centesis Catheter Needle	17g	12cm	n/a	n/a	Clear Hub	n/a	1 ea
NDL-009-06	17g x 17cm, Centesis Catheter Needle	17g	17cm	n/a	n/a	Clear Hub	n/a	1 ea
NDL-009-01	19g x 11cm, Centesis Catheter Needle	19g	11cm	n/a	n/a	Clear Hub	n/a	1 ea
NDL-009-03	19g x 17cm, Centesis Catheter Needle	19g	17cm	n/a	n/a	Clear Hub	n/a	1 ea
NDL-009-02	21g x 11cm, Centesis Catheter Needle	21g	11cm	n/a	n/a	Clear Hub	n/a	1 ea

Chiba Needles

SKU	Product Description	Gauge	Length	Non Echo/Echo	Bevel	Hub Color	Additional Features	UoM
NDL-013-00	18g x 9cm, Chiba Needle	18g	9cm		n/a	n/a	n/a	1 ea
NDL-011-00	22g x 15cm, Chiba Needle	22g	15cm	Echo	n/a	n/a	Depth Markers	1 ea
NDL-011-01	22g x 20cm, Chiba Needle	22g	20cm	Echo	n/a	n/a	Depth Markers	1 ea

Custom sizes and modifications available | Can be sold in either BNS and/or sterile private label



GALT™

OEM Division

A Theragenics Company

To order:

phone 1-972-271-5177 or 1-800-639-2800

email: galtorders@galtmedical.com

Galt Medical Corp. 2220 Merritt Drive, Garland, TX 75041

phone: 1-972-271-5177 email: info@galtmedical.com

www.galtmedical.com