

Indicated for Neonatal  
and Pediatric Use


# MicroSlide®

Pediatric MST Products

For smaller patient access, keep it short and sweet.

Galt Medical Corp. introduces the MicroSlide® line of vascular access products for smaller patients and smaller vessels. Maintaining the features that are preferred in our current offerings, this gives the clinician an added advantage to size the product to the patient. Patients benefit from being accessed by the appropriately sized products. It is a win-win solution, and now available from Galt Medical.

- 2Fr and 3Fr tearaway introducers
- Available individually or in complete MST Kit
- .008", .010" and .014" diameter guidewires
- 24ga x 2.5cm and 24ga x 4cm needles



(Actual Size)



# MicroSlide®

## Pediatric MST Products

### Growth Stages

The MicroSlide® Pediatric MST product line is the first to offer both a 2Fr and 3Fr MST offering. Galt continues to offer the flexibility of a full MST Kit or sterile, separate components. Our MST Tearaway line offering continues to grow and now includes 3cm, 5cm, 7cm and 10cm lengths. Our popular .008" and .010" guidewire configurations are now offered in 20cm, 40cm and 60cm lengths. The addition of the .014" offering continues the commitment of solution-based products for clinicians and patients.



(Actual Size)

#### Tearaway Introducer Sets

Set Contents: PTFE Sheath and Dilator only

Catalog Number	French Size	Length	Accepts Guidewire	Units/Box
INT-101-02	2F	3cm	.008"	10
INT-101-14	2F	3cm	.014"	10
INT-101-15	2F	7cm	.014"	10
INT-101-03	3F	3cm	.018"	10



#### Tearaway Introducer Kits

Catalog Number	French Size	Length	NT Guidewire	Echo Needle	Units/Box
KIT-051-00	2F	3cm	.008" x 20cm	24ga x 2.5cm	5
KIT-051-01	2F	3cm	.008" x 40cm	24ga x 2.5cm	5
KIT-051-02	2F	3cm	.008" x 20cm	24ga x 4cm	5
KIT-051-03	2F	3cm	.008" x 40cm	24ga x 4cm	5
KIT-051-04	3F	3cm	.010" x 20cm	24ga x 2.5cm	5
KIT-051-05	3F	3cm	.010" x 40cm	24ga x 2.5cm	5
KIT-051-06	3F	3cm	.010" x 20cm	24ga x 4cm	5
KIT-051-11	3F	3cm	.010" x 40cm	24ga x 4cm	5

#### Guidewires

Nitinol mandrel-tapered both ends

Catalog Number	Size	Length	Units/Box
SGW-051-06	.008"	20cm	5
SGW-051-08	.008"	40cm	5
SGW-051-03	.008"	60cm	10
SGW-051-07	.010"	20cm	5
SGW-051-09	.010"	40cm	5
SGW-051-04	.010"	60cm	10

#### NT/Gold Wire

Nitinol mandrel-gold tip guidewires

Catalog Number	Size	Length	Units/Box
SGW-056-62	.014"	20cm	5
SGW-056-60	.014"	40cm	5



(Actual Size)

#### Needles

Catalog Number	Size	Length	Tip	Units/Box
NDL-100-29	24ga	2.5cm	echogenic	25
NDL-100-30	24ga	4cm	echogenic	25

