




## Valved Tearaway Introducers



Hemostasis Valve



Smooth Transition

### Reliable. Consistent. Low Insertion Force.

The GaltVTI valved tearaway introducer sets feature an integrated, low insertion force hemostasis valve, balancing decreased risk of blood loss and air embolism with ease of catheter insertion.

#### Features:

- Ergonomic design with smooth transition
- Easy, secure dilator locking mechanism
- Consistent peel
- Available in 5Fr to 16Fr



For more information about GaltVTI Valved Tearaway Introducers, call (800) 639-2800 or visit [galtmedical.com](http://galtmedical.com).

## GaltVTI® Valved Tearaway Introducer Sets

Model Number	Size (Fr)	Sheath Length (cm)	Dilator Length (cm)	Accepts Guidewire Diameter	Box Quantity
VTI-007-05	5	13	17	.038"	10
VTI-007-06	6	13	17	.038"	10
VTI-007-07	7	13	17	.038"	10
VTI-007-08	8	13	17	.038"	10
VTI-007-09	9	13	17	.038"	10
VTI-007-10	10	13	17	.038"	10
VTI-007-11	11	13	17	.038"	10
VTI-007-12	12	13	17	.038"	10
VTI-007-13	13	13	17	.038"	10
VTI-007-14	14	13	17	.038"	10
VTI-007-15	15	13	17	.038"	10
VTI-007-16	16	13	17	.038"	10

