

ELITE|HV T

TRANSITIONLESS HEMOSTASIS
VALVE INTRODUCER



**UP TO 40% LESS ENTRY/INSERTION
FORCE VS COMPARABLE SHEATH.***

*1:1 simulated use test data on file at Galt Medical Corp.

The Elite HV T .018" combination micro-introducer and valve introducer sheath delivers transitionless insertion force and over 40% improvement to current users. This compact all-in-one sheath system reduces steps, saves time, and minimizes product inventory.

✉ info@galtmedical.com

☎ Phone: (800) 639-2800 - Fax: (972) 271-4706



Galt Medical Corp. - Garland, Texas




www.galtmedical.com



ELITE|HV[®] 7

TRANSITIONLESS HEMOSTASIS VALVE INTRODUCER



✓ Available in set or kit configurations
(set = dilator/sheath assembly only).

✓ Increased 4Fr and 5Fr ID to accommodate
oversized catheter size variations.

✓ Kits include multiple needle lengths to
size correctly according to patient.

	Catalog No.	Sheath Size	Sheath Length	Needle Length	R/O Band	Accepts Guidewire Diameter	Guidewire
Set*	CLI-109-45	4F	7cm	N/A	No	.018"	N/A
Set*	CLI-109-55	5F	7cm	N/A	No	.018"	N/A
Set*	CLI-109-06	6F	7cm	N/A	No	.018"	N/A
Set*	CLI-109-07	7F	7cm	N/A	No	.018"	N/A
Kit	CLI-122-45	4F	7cm	2.5cm/4cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-122-55	5F	7cm	2.5cm/4cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-122-06	6F	7cm	2.5cm/4cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-122-07	7F	7cm	2.5cm/4cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-124-45	4F	11cm	4cm/7cm	No	.018"	Stainless Mandrel & Tip
Kit	CLI-124-55	5F	11cm	4cm/7cm	No	.018"	Stainless Mandrel & Tip
Kit	CLI-124-06	6F	11cm	4cm/7cm	No	.018"	Stainless Mandrel & Tip
Kit	CLI-124-07	7F	11cm	4cm/7cm	No	.018"	Stainless Mandrel & Tip
Kit	CLI-124-14	4Fr	11cm	4cm/7cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-124-15	5Fr	11cm	4cm/7cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-124-16	6Fr	11cm	4cm/7cm	No	.018"	Nitinol Mandrel/Gold Tip
Kit	CLI-124-17	7F	11cm	4cm/7cm	No	.018"	Nitinol Mandrel/Gold Tip

Note: 10 Units/Box

*Set includes sheath/dilator assembly only.

Data on file at Galt Medical

✉ info@galtmedical.com

☎ Phone: (800) 639-2800 - Fax: (972) 271-4706



Galt Medical Corp. - Garland, Texas



www.galtmedical.com

