

Guide Wire Torque Device



One-handed control

The quick press and release handling system offers an easy and efficient alternative to the current torque device options.

Two options available:

- The yellow torque device is used with guidewires with diameters between .010" and .020"
- The orange torque device is used with guidewires with diameters between .025" and .040"

Instructions for use:

1. Press open end firmly and insert over guidewire.
2. Position torque device and release.
3. Press finger grips to adjust or remove.

Ordering information:

TOR-000-10 Accepts guidewire diameter .010"-.020" Yellow

TOR-000-25 Accepts guidewire diameter .025"-.040" Orange

MKT-100-98 Orig.

