


Hemostasis Valve Adapter



Hemostasis Valve Adapter

The attachable hemostasis valve adapter is the answer to minimizing blood loss during the introduction of catheters or devices into the vasculature. For those times in procedures when a device needs to “park” for a period of time, the Valve Adapter is your go-to device.

The Valve Adapter consists of a luer fit hub, a hemostasis valve, and side arm tubing with a stopcock.

Features:

- Ability to pass up to an 8Fr device
- Ability to seal on mandrel wires
- Available bulk non-sterile and sterile private label
- Pressure rated up to 30 PSI

Ordering Information:

ACL-000-00 Hemostasis Valve Adapter 10/box



For more information about the Valve Adapter call (800) 639-2800 or visit galtmedical.com.

MKT-100-86 Orig.